
Rev From	Rev To	ECO #	Date	Appr.
0599	0803	2021-03		

Notes

- BD Cat. No. 276709
- Blank (Sheet) Size : Length: NA Width: NA
 Number of Pages: 8 Number of Sheets: 1
 Page Size: Length 16" Width 3.625" Final Folded Size: NA
- Style (see illustrations below): # NA

- See Specification Control No. NA for Material Information
- Ink Colors: Printed two sides Yes No
 No. of Colors: 1 Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS CONTROLLED BY BD CARIBE, LTD.

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA	
Proofer	Date			
Checked By	Date			
Part Number:	0210745	Category and Description Package Insert MacroVue Control	Sheet: 1 of 9 Scale:	A

BD Macro-Vue™ RPR Card Test Control Cards

0210745
2003/08

English: pages 1 - 2 Italiano: pagine 5 - 6
Français : pages 2 - 3 Español: páginas 6 - 7
Deutsch: Seiten 3 - 5

See symbol glossary at end of insert. / Se symbolglossaret i slutningen af indlægssedlen. / Voir le glossaire des symboles à la fin de la notice. / Siehe Symbol-Erklärungen am Ende der Packungsbeilage. / Δείτε το γλωσσάριο των συμβόλων στο τέλος του ένθετου. / Vedere il glossario dei simboli alla fine del foglio illustrativo. / Consulte o glossário de símbolos no fim do folheto informativo. / Consulte el glosario de símbolos al final del prospecto. / Se symbolförteckningen vid slutet av bipacksedeln.

Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Contacte o seu representante local da BD para obter instruções. / Kontakta lokal Becton Dickinson-representant för anvisningar.

INTENDED USE

The **Macro-Vue™** RPR Card Test Control Cards are dehydrated control specimens of predetermined reactivity for quality-control testing of antigen before use in performing the **Macro-Vue** RPR Card Tests for the serologic detection of syphilis.

SUMMARY AND EXPLANATION

Macro-Vue RPR Card Test Control Cards were designed for routine, daily use as check-test specimens to determine the acceptable reactivity of RPR Card antigen suspension. RPR Card antigen is a cardiolipin type antigen, which detects "reagin", an antibody-like substance present in serums from syphilitic persons, and occasionally in serums of persons with other acute and chronic conditions.¹⁻⁴

PRINCIPLES OF THE PROCEDURE

The RPR 18 mm Circle Control Cards have labelled test areas containing dehydrated specimens with designated patterns of graded reactivity. Upon rehydration and proper execution of the Control Card procedure, the RPR Card antigen suspension should reproduce the established reactivity pattern (see "Procedure") to confirm optimal reactivity of the antigen.

REAGENTS

The RPR 18 mm Circle Card Test Control Card consists of three circles designated as: Reactive, Reactive minimal-to-moderate and Nonreactive. The labelled areas on this card contain graded reactivity specimens. The specimen material is human plasma purchased from licensed blood banks as Reactive or Nonreactive for syphilis.

Warnings and Precautions:

For *in vitro* Diagnostic Use.

Before proceeding, review "Procedure" and "Warnings and Precautions" in the package insert supplied in all **Macro-Vue** RPR Card Test Kits.

The plasma used on the control cards is derived from human blood tested by an FDA (U.S. Food and Drug Administration)-approved method for the presence of the antibody to HIV (human immunodeficiency virus) and HBsAg (hepatitis B surface antigen) and found to be nonreactive.

WARNING: Because no test method can offer complete assurance that HIV, hepatitis B virus, or other infectious agents are absent, SPECIMENS AND THESE REAGENTS SHOULD BE HANDLED AS THOUGH CAPABLE OF TRANSMITTING AN INFECTIOUS DISEASE. "Standard Precautions"⁵⁻⁸ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Controls and RPR Card antigen suspension should be at room temperature (23 to 29°C) when used.

If the established reactivity pattern is not obtained, delay routine testing until the optimal reactivity is obtained by rechecking antigen suspension, test procedures, room temperature and equipment (including use of a rotator with humidifying cover). Non-acceptable results invalidate tests on individual specimens.

Once the dehydrated specimens on the control card are reconstituted, the specimen solutions should be tested immediately using the RPR Card Test procedure. Delay in testing may allow evaporation to affect the reactivity pattern.

Prolonged storage at cool room temperatures (21 to 25°C) and at elevated temperatures (above 27°C) causes deterioration of the specimens, showing a decrease in reactivity and difficulty in reconstituting the card. Once the sealed envelope is opened, the card should be used immediately.

Reconstitution of Control Card: Each test area of the control card should be reconstituted with 0.05 mL of *deionized/distilled water*.

Storage Instructions: In the unopened foil-sealed envelopes, the control cards retain an acceptable reactivity pattern for the assigned expiration date when stored under refrigeration at 2 to 8°C, which is the recommended storage temperature.

PROCEDURE AND RESULTS

Materials Provided: RPR Card Test Control Cards

Materials Required But Not Provided:

- One of the following **Macro-Vue** RPR Card Test Kits*:
18 mm Circle Test—Kit No. 115, 110, 104, 510 or 532 (qualitative) or Kit No. 112 (quantitative)

All kits contain sufficient coverage of materials for performance of daily control card tests.

- A micropipettor which dispenses 0.05 mL (50 µL) and 0.03 mL (30 µL).
- A rotator, 100 ± 2 rpm, 2 cm diameter, automatic timer, friction drive, with a humidifying cover containing a moistened blotter or sponge.

*For availability of kits, refer to the package insert supplied with all **Macro-Vue** RPR Card Test Kits.

Test Procedure:

- Place 0.05 mL (50 µL) deionized/distilled water onto the three circles, using an automatic micropipettor (or 0.05 mL capillary with rubber bulb attached). *Do not use **Dispenstirs™** device to reconstitute a Control Card.*
- Using the stirring end of an 0.05 mL **Dispenstirs** device (or broad end of a stirrer) for each circle, mix until dehydrated control specimen is dissolved. Spread specimen to the inside of the circle. *A new **Dispenstirs** device or stirrer must be used for each circle!*
- Shake the antigen bottle before use. Holding in a vertical position, dispense one or two drops in the cap to make sure needle passage is clear. Place one "free-falling" drop of antigen (1/60 mL antigen – [20 G yellow hub needle]) onto each test area. *Do not restir;* mixing of antigen and specimen is accomplished during rotation. Pick up pre-dropped antigen from the bottle cap.
- Rotate for eight min, on a mechanical rotator at 100 ± 2 rpm, using a moistened humidifying cover. Following rotation, to help differentiate Nonreactive from Reactive minimal-to-moderate results, a brief rotating and tilting of the card by hand (three or four to-and-fro motions) must be made.
- Immediately* read the card macroscopically in the "wet" state under a high intensity incandescent lamp or strong daylight.

The Reactive control should show characteristic strong clumping; the Nonreactive control should show the smooth, grayish appearance of unclumped particles. The Reactive minimal-to-moderate control should show minimal-to-moderate clumping. See the Reading Guide in the package insert provided in all RPR Card Test Kits.

User Quality Control:

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

LIMITATIONS OF THE PROCEDURE

If the specimens do not dissolve satisfactorily when reconstituted with deionized/distilled water, the card should not be used. Delay in performance of the test procedure following reconstitution, or performing the test procedure in an environment of low humidity (particularly without a humidifying cover on the rotator) may cause an undesirable reactivity pattern (a granular appearance in the Nonreactive specimen). Since the dehydrated specimens are quantitatively controlled for reactivity, the proper execution of the Control Card and RPR Card Test procedures are essential for the determination of the expected reactivity pattern.

The intensity of clumping observed should not be used as a reading standard or guide in interpreting reactions observed with unknown specimens. The Reading Guide contained in the package insert provided with all RPR Card Test Kits is intended for this purpose.

EXPECTED VALUES AND PERFORMANCE CHARACTERISTICS

Each lot of 18 mm Circle Control Cards is tested and compared for the designated pattern of reactivity against reference antigen suspensions and meets the U.S. Centers for Disease Control and Prevention (CDC) product specifications.

The expected patterns of reactivity on the control cards following the performance of RPR Card Test procedures are shown in the Reading Guides contained in the package insert provided in all **Macro-Vue** RPR Card Test Kits.

Reactive Circle – Similar to undiluted Reactive titer spot shown under the RPR Card quantitative test.

Reactive minimal-to-moderate Circle – Similar to Reactive minimal spot of RPR Card qualitative test and to Reactive (1:4) spot of RPR Card quantitative test.

Nonreactive Circle – Similar to RPR Card Nonreactive spot.

AVAILABILITY**Cat. No. Description**

276709 **Macro-Vue™** RPR Card Test Control Cards,
Box of ten individual foil-sealed cards:
18 mm Circle (Reactive, Nonreactive and Reactive Minimal
to Moderate circles), Box of 10.

REFERENCES

- Creighton, E.T. 1990. Rapid plasma reagin (RPR) 18-mm circle card test, p. 99-108. *In* S.A. Larsen, E.F. Hunter, and S.J. Kraus (ed.), A manual of tests for syphilis, 8th ed. American Public Health Association, Washington, D.C.
- Portnoy, J., J.H. Brewer, and A. Harris. 1962. Rapid plasma reagin card test for syphilis and other treponematoses. *Public Health Rep.* 77:645-652.
- Portnoy, J. 1963. Modifications of the rapid plasma reagin (RPR) card test for syphilis, for use in large scale testing. *Am. J. Clin. Pathol.* 40:473-479.
- Falcone, V.H., G.W. Stout, and M.B. Moore, Jr. 1964. Evaluation of rapid plasma reagin (circle) card test. *Public Health Rep.* 79:491-495.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262, 17/10/2000, p. 0021-0045.*

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BENEX Limited