

For more information contact the Infolab Distribution Center nearest you:

Clarksdale, MS	1-800-647-8222
Fort Oglethorpe, GA	1-800-241-8953
Greensboro, NC	1-800-782-9700
Round Rock, TX	1-800- 322-3882
Ocala, FL	1-800- 734-7584
Zionsville, IN	1-866- 232-7096



**Ac-T Diff 2** delivers performance, value and safety for physician office labs and smaller laboratories. It's an affordable, easy-to-use system that offers a small sample volume in both open and closed analysis modes, a wide operating range and multiple output formats.



**Beckman Coulter AcT Diff II Hematology Analyzer**

**Specifications:**

- **Closed Vial Sampling**
- **Only 18 µL Sample Size**
- **Compact**
- **Zero routine maintenance**

**Imprecision Specifications**

Parameter	Level	Units	CV%	SD
WBC	6.0 - 15.0	x10 <sup>3</sup> cells/µ	<3.0	
RBC	3.0 - 6.0	x10 <sup>6</sup> cells/µ	<3.0	
Hgb	12.0 - 18.0	g/dL	<2.0	
MCV	80.0-100.0	fL	<3.0	
Plt	200 - 500	x10 <sup>3</sup> cells/µ	<7.0	
MPV	5.0 - 20.0	fL	<3.0	
RDW	12.0 - 15.0	%	<3.0	
LY	20 - 50	%		<1.5
MO	2.0 - 10.0	%		<1.5
GR	30.0 - 70.0	%		<3.0

**Operating Range**

Parameter	Range	Units
WBC	0.0 150	x10 <sup>3</sup> cells/µ
RBC	0.00 8.00	x10 <sup>6</sup> cells/µ
Hgb	00.0 - 30.0	g/dL
MCV	50.0 - 130.0	fL
Plt	000 - 3000	x10 <sup>3</sup> cells/µ
MPV	5.0 - 20.0	fL
LY	0 - 100	%