



# HEMOCHRON®

Whole Blood Coagulation Systems

*Signature+*

## FEATURES:

- Performs rapid and accurate ACT+, ACT-LR, APTT, and PT tests within minutes
- Customize the Signature+ instrument with programming options including:
  - two level QC lockout
  - preprogrammed valid operator IDs and notes
- Store up to 400 QC results and 400 patient results.
- Offers liquid and electronic quality controls



## AVAILABLE TESTS:

**ACT-LR** - Rapidly monitors low to moderate levels of heparin (0-2.5 units/ml), reports results as Celite® equivalent values.

**ACT+** - Excels at monitoring moderate to high levels of heparin (1-6 units/ml). Unaffected by aprotinin, hypothermia, and hemodilution. Results are displayed as Celite- equivalent values.

**APTT** - Displays sensitivity to low-dose heparin anticoagulation up to 1.5 units/ml. Reports plasma-equivalent values up to 400 seconds.

**PT** - Provides the flexibility to utilize venous or fingerstick samples. Reports results as International Normalized Ratio (INR) and as plasma-equivalent values.

**Citrate PT** - For PT testing of a citrated whole blood sample. Reports results as International Normalized Ratio (INR) and as plasma-equivalent values.

**Citrate APTT** - For APTT testing of a citrated whole blood sample. Reports plasma-equivalent values up to 400 seconds.





## CONFIGURATION MANAGER

Customize each Signature+ to your individual needs with numerous available options.

### Operator ID and Patient ID

- Program up to 100 valid operator IDs.
- Lockout capabilities allow access to specific system options.
- ID can be required for patient testing and/or quality control testing.
- Reuse IDs for faster tagging of multiple tests performed on one patient, which is ideal for reducing the number of steps to run tests during long procedures.

### Quality Control Lockout

- Ensures timely QC testing and CLIA compliance for both liquid and electronic quality control.
- Lockout by time lapse and/or QC test failure.
- “911” option allows and tracks override of lockout in emergency situations.
- Number of 911s controlled by supervisor.

### Security

- Password protect instrument by operator.
- Configuration Manager password protected by supervisor.
- Limited access to set-up and control options reduces accidental setting changes.
- Prohibit users from changing date and time or purging database.

### Customized Notes

- Attach up to two preloaded notes to patient results with one touch.
- Print with results for hard copy explanations and comments.
- Scroll selected notes on instrument.

### Compatible with Windows® 95/98/2000/NT/XP



Clarksdale, MS  
Toll Free: 1-800-647-8222  
Phone: (662) 627-2283  
Fax: (662) 627-1913  
Email:  
[infolabms@infolabinc.com](mailto:infolabms@infolabinc.com)

Round Rock, TX  
Toll Free: 1-800- 322-3882  
Phone: (512) 255-7994  
Fax: (512) 244-9659  
Email:  
[infolabtx@infolabinc.com](mailto:infolabtx@infolabinc.com)

Fort Oglethorpe, GA  
Toll Free: 1-800-241-8953  
Phone: (706) 861-5046  
Fax: (706) 866-5234  
Email:  
[infolabga@infolabinc.com](mailto:infolabga@infolabinc.com)

Ocala, FL  
Toll Free: 1-800- 734-7584  
Phone: (352) 237-9654  
Fax: (352) 237-9733  
Email:  
[infolabfl@infolabinc.com](mailto:infolabfl@infolabinc.com)

Greensboro, NC  
Toll Free: 1-800-782-9700  
Phone: (336) 272-2874  
Fax: (336) 274-9473  
Email:  
[infolabnc@infolabinc.com](mailto:infolabnc@infolabinc.com)

Zionsville, IN  
Toll Free: (866) 232-7096  
Phone: (317) 873-4400  
Fax: (317) 873-4443  
Email:  
[infolabin@infolabinc.com](mailto:infolabin@infolabinc.com)