



COULTER® LH 750 Hematology Analyzer




Performance and Results that Exceed Your Expectations

The LH 750 analyzer is equipped with superior technologies, NRBC enumeration and random access capabilities. These advancements help high-volume labs eliminate manual pre-sorting of samples and speed results to clinicians. With this system, you can improve overall efficiency, while maintaining the highest degree of testing accuracy.



Product Features

COULTER® LH 750

WBC LINEARITY



PLATELET LINEARITY




Improve Cell Enumeration Accuracy With AccuCount Technology

- Improve productivity with WBC linearity that extends from zero to 400,000 and platelet linearity that extends from zero to 3,000,000
- Increase clinician confidence at clinical decision thresholds without additional costly testing protocols
- Enhance productivity with a reduction in false-positive results and confirmatory tests

Simplify and Automate NRBC Enumeration

- Eliminate the use of additional reagents since NRBCs are automatically enumerated with every CBC/diff

- Identify and enumerate NRBCs without manual intervention
- Improve productivity by reducing the need for technologist intervention and sample re-runs
- Obtain precise and accurate WBC counts with automatic interference correction



Streamline Data Management with Improved Decision Support

- Standardize laboratory processes across all shifts with user-defined decision rules
- Customize the system to meet your lab's needs and specific patient populations with flexible differential flagging



Receive Unparalleled Service and Support

- Increase uptime and productivity with on-board diagnostics and troubleshooting assistance, or with PROService Remote Diagnostics support*
- Receive top-quality customer support through our telephone, on-line or on-site services

More Info

Increase Testing Confidence with AccuCount Technology

AccuCount Technology is our revolutionary technique that combines the power of the Coulter Principle with our most advanced algorithms to deliver the best overall precision, linearity and accuracy in the industry. With greatly reduced false positive flagging, AccuCount Technology dramatically streamlines your testing process. It eliminates unnecessary operator intervention such as manual WBC correction for making blood smears, sample dilution preparations and repeat testing.

* PROService is only available in the U.S.

Note: Not all products are available in all countries.



For More Information Contact:

The Infolab Distribution Center Nearest You

Clarksdale, MS**Toll Free: 1-800-647-8222****Phone: (662) 627-2283****Fax: (662) 627-1913****Email:****infolabms@infolabinc.com****Fort Oglethorpe, GA****Toll Free: 1-800-241-8953****Phone: (706) 861-5046****Fax: (706) 866-5234****Email:****infolabga@infolabinc.com****Greensboro, NC****Toll Free: 1-800-782-9700****Phone: (336) 272-2874****Fax: (336) 274-9473****Email:****infolabnc@infolabinc.com****Round Rock, TX****Toll Free: 1-800- 322-3882****Phone: (512) 255-7994****Fax: (512) 244-9659****Email:****infolabtx@infolabinc.com****Ocala, FL****Toll Free: 1-800- 734-7584****Phone: (352) 237-9654****Fax: (352) 237-9733****Email:****infolabfl@infolabinc.com****Zionsville, IN****Toll Free: (866) 232-7096****Phone: (317) 873-4400****Fax: (317) 873-4443****Email:****infolabnc@infolabinc.com**