

VITROS® 950 Chemistry System

Reliably accurate results
for high-volume testing.




With the VITROS 950, you're guaranteed that the result you get is the right result. Anyplace, anytime. The system requires no reagent preparation, and with no plumbing and minimal maintenance, there's virtually no downtime. Just increased productivity and reduced operating costs, so you spend less time re-running tests and have less need for a back-up system. And the system's high throughput, large capacity and expanding menu go a long way to meeting the demands of today's high-volume workload. So take a closer look and you'll see that the VITROS 950 is the best choice for Value and Innovation from a company you can Trust.

Distributed by Infolab, Inc. 1-800-782-9700

VITROS® Do More. For Life.

Superior results. High-volume reliability. Cost-effective operation.

The system of choice for consistently reliable, high-quality results.



- Proven MicroSlide™ technology sets the standard for superior assay precision
- Minimal interference from hemolysis, lipemia and bilirubin

- Fewer repeats, redraws and dilutions than with wet chemistry systems
- Backed by the company top rated in service and support for over 17 years

The best checks and balances on any chemistry system.

- Short sample detection/clot detection/bubble detection—ensure accuracy
- Unique Smart metering™—ensures dependable testing integrity
- Disposable tips—ensure no carryover
- Continuous process verification—ensures the right result the first time

The power to meet the demands of a high-volume workload.

- Throughput of up to 900 results/hour
- 40+ on-board chemistries
- Easily upgraded to connect to any automation sample transport system

Easy to operate—just plug in and go.

- No deionized water, no plumbing, no drains
- Load bar-coded slide cartridges, of any lot number, in one easy step—with no reagent preparation
- Accepts any size primary/secondary containers or cups
- Small sample size requirements
- User-friendly help text is available at all times
- Biosafe operations with no infectious waste

Reduce your costs. Streamline your workflow.

- The highest reportable result efficiency of any chemistry analyzer—90-95%
- Unit-dose test slides provide predictable operating costs
- QC once daily/up to 6-month calibration stability for most assays
- Minimal maintenance requirements
- Consistent operator interface with other VITROS systems—operate one, operate all

For more information, contact your Ortho-Clinical Diagnostics Representative or call 1-800-828-6316, prompt #1.

Distributed by Infolab, Inc. 1-800-782-9700



Ortho-Clinical Diagnostics

a Johnson & Johnson company