



## IMMULITE® 1000

Key features of the IMMULITE 1000 Immunoassay System:



### Reliable Performance

Maintains the strong features of IMMULITE: extensive menu, ease of use, reliable and robust hardware, and excellent assay performance.

### Comprehensive Menu

The breadth of menu on the IMMULITE 1000 provides an opportunity for true consolidation of all tests for a given patient on one analyzer, resulting in cost savings and laboratory efficiency.

### IMMULITE Allergy Testing

IMMULITE 1000 assays distinguish atopic from nonatopic individuals and identify the class of allergen: grass, weed, tree, animal, mold or food. The fast turnaround time produces latex-sensitivity results in just 72 minutes.

### True Random Access Capability

Any test or combination of tests can be run on the IMMULITE 1000 at any time. Samples can be continuously processed for improved turnaround times. STAT samples can also be added easily whenever required so that the processing time is kept to a minimum.

### Consolidated Footprint

Computer, keyboard, monitor, and built-in probe wash and water containers integrated into the system, for reduced system footprint.

### Windows®-Based Software

State-of-the art, icon-driven Microsoft® Windows software designed for intuitive, minimal user interaction; operation via touch screen, trackball or keyboard. Excellent QC package with online QC tracking and trending (Levey–Jennings plot). Flexible and customizable patient reports. Multilingual support for six languages.

### Onboard Dilutions

Ability to perform onboard dilutions for infectious disease and other selected prediluted assays, as well as HCG.

### After-Sales Service

DPC is proud to supply all our customers with after-sales service that they can rely on. We are committed through our Quality Standards to support our customers in every way possible and we continually strive to improve the service and support we offer.

### Positive Sample ID

Modified sample cup holder that allows for running a barcoded secondary microtube on board, for positive secondary sample identification.

### Remote Diagnostics

This feature provides customers with enhanced support and troubleshooting assistance that can reduce the downtime of the system.

### Easy-to-Use Design

The benchtop design can easily be incorporated into any laboratory area. The software has been developed to maximize the flexibility and productivity of the laboratory.



---

## For More Information Contact:

### The Infolab Distribution Center Nearest You

**Clarksdale, MS****Toll Free: 1-800-647-8222****Phone: (662) 627-2283****Fax: (662) 627-1913****Email:****[infolabms@infolabinc.com](mailto:infolabms@infolabinc.com)****Fort Oglethorpe, GA****Toll Free: 1-800-241-8953****Phone: (706) 861-5046****Fax: (706) 866-5234****Email:****[infolabga@infolabinc.com](mailto:infolabga@infolabinc.com)****Greensboro, NC****Toll Free: 1-800-782-9700****Phone: (336) 272-2874****Fax: (336) 274-9473****Email:****[infolabnc@infolabinc.com](mailto:infolabnc@infolabinc.com)****Round Rock, TX****Toll Free: 1-800- 322-3882****Phone: (512) 255-7994****Fax: (512) 244-9659****Email:****[infolabtx@infolabinc.com](mailto:infolabtx@infolabinc.com)****Ocala, FL****Toll Free: 1-800- 734-7584****Phone: (352) 237-9654****Fax: (352) 237-9733****Email:****[infolabfl@infolabinc.com](mailto:infolabfl@infolabinc.com)****Zionsville, IN****Toll Free: (866) 232-7096****Phone: (317) 873-4400****Fax: (317) 873-4443****Email:****[infolabnc@infolabinc.com](mailto:infolabnc@infolabinc.com)**